

Accuracy Range of Thermometer (+/-) _____ °C
*** Record temperatures at least twice per day ***
Fridge: _____
** Must write in numbers to one decimal place – DO NOT Use (X) or (✓)*

Day of the Month	1		2		3		4		5		6		7		8		9		10		11		12		13		14		15	
Exact Time																														
	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm	am	pm
Refrigerator Temperature °C	≥ 10	Take immediate ACTION if ANY temperature reading is in a shaded section (Current, Min or Max) * Consider the <i>accuracy range of thermometer</i> when plotting temperatures for possible Cold Chain Incidents (CCIs) * Ex → 7.8°C + 0.5°C (+/-) range = Possible temp 8.3°C (CCI)																												
	9																													
	8																													
	7																													
	6																													
	5																													
	4																													
	3																													
	2																													
	1																													
	0																													
≤ -1	Take immediate ACTION if ANY temperature reading is in a shaded section * Consider the <i>accuracy range of thermometer</i> when plotting temperatures for possible Cold Chain Incidents (CCIs) *																													
Current °C																														
Min Temp °C																														
Max Temp °C																														
Room °C																														
Dial / Temp Setting																														
Initials																														

Write Current, Min and Max Temps with 1 decimal place

If vaccine has been exposed ANY TIME to temperatures outside of +2°C to +8°C (after factoring in the thermometer accuracy range) follow the Cold Chain Incident guidelines and recommendations

- Current Temp** → The current “now” temperature of the refrigerator – This shows a small window of time
- Min / Max Temp Range** → The highest and lowest temperature reached *in the total time since the thermometer was last reset* - **The most rich and important temperature data**
- Dial/Temp Setting** → The programed temperature or dial setting on the fridge - **The recommended programed setting should allow the fridge’s resting temperature to sit between +4.0°C to +5.0°C**

Refrigerator Temperature Form (°C)

Month _____ Year _____

Accuracy Range of Thermometer (+/-) _____ °C

* Record temperatures at least twice per day *

Fridge: _____

 * **Must write in numbers to one decimal place** – DO NOT Use (X) or (✓)

Day of the Month		16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Exact Time		am pm		am pm		am pm		am pm		am pm		am pm		am pm		am pm	
Refrigerator Temperature °C	≥ 10	Take immediate ACTION if ANY temperature reading is in a shaded section (Current, Min or Max) * Consider the <i>accuracy range of thermometer</i> when plotting temperatures for possible Cold Chain Incidents (CCIs) * Ex → 7.8°C + 0.5°C (+/-) range = Possible temp 8.3°C (CCI)															
	9																
	8																
	7																
	6																
	5																
	4																
	3																
	2																
	1																
0	Take immediate ACTION if ANY temperature reading is in a shaded section * Consider the <i>accuracy range of thermometer</i> when plotting temperatures for possible Cold Chain Incidents (CCIs) *																
≤ -1																	
Current °C																	
Min Temp °C																	
Max Temp °C																	
Room °C																	
Dial / Temp Setting																	
Initials																	

Write Current, Min and Max Temps with 1 decimal place

 If vaccine has been exposed ANY TIME to temperatures outside of +2°C to +8°C (**after factoring in the thermometer accuracy range**) follow the Cold Chain Incident guidelines and recommendations

Current Temp → The current “now” temperature of the refrigerator – This shows a small window of time

Min / Max Temp Range → The highest and lowest temperature reached *in the total time since the thermometer was last reset* - **The most rich and important temperature data**
Dial/Temp Setting → The programmed temperature or dial setting on the fridge - **The recommended programed setting should allow the fridge’s resting temperature to sit between +4.0°C to +5.0°C**